

a true bathing experience

Comparison Chart

	Endres Bellavita	Mountway Aquila	Minivator 303	Invacare Orca	Aqua Joy Premier	Homecraft Sonaris
Bathlift Weight	9.3kg	9.6kg	13.8kg	12.1kg	12kg	12kg
Minimum Seat Height (Lowered)	6cm	7.5cm	7.5cm	6cm	8cm	7.5cm
Maximum Seat Height	48cm	44cm	41cm	42cm	46cm	45cm
Angle of Recline	50°	35°	45°	40°	40°	40°
Maximum Weight Capacity	140kg	140kg	135kg	135kg	170kg	140kg
Hygienic Cutaway	Yes	No	No	No	No	No
NEW Lithium Battery Technology	Yes	No	No	No	No	No
Anti-Bacterial treated covers as standard	Yes	No	Yes	Yes	Yes	Yes
Bathlift Guarantee	Lifetime Guarantee on Frame & Motor	5 Yrs on Frame	2 Yrs	3 Yrs on Frame & Parts	3 Yrs on Frame & Motor	3 Yrs on Frame
	2 Yrs on Handset Charger & Battery	2 Yrs on Battery & Parts			2 Yrs on Handset & Battery	2 Yrs on Handset & Charger



Telephone: **01422 314488**
 Email: **enquiries@drivemedical.co.uk**
www.drivemedical.co.uk

drive ^{CE}
 MEDICAL DESIGN & MANUFACTURING
www.drivemedical.co.uk