



OUR EXCELLENT RANGE OF POWERED MOBILITY PRODUCTS CONTINUES TO DEVELOP!



Geo Micro Code: S3500

LENGTH 85cm (34")	HEIGHT 85cm (34")	WIDTH 56cm (22")
SEAT WIDTH 45cm (18")	MAX CAPACITY 115kg (18St)	WEIGHT 46kg (101lb)
MAXIMUM GRADIENT 6°	MAXIMUM RANGE 6 Miles (10km)	

Sunfire Plus GT Code: PCMM07

LENGTH 104cm (41")	HEIGHT 114cm (45")	WIDTH 66cm (26")
SEAT WIDTH 45cm (18")	MAX CAPACITY 136kg (21St)	WEIGHT 96kg (201lb)
MAXIMUM GRADIENT 8°	MAXIMUM RANGE 20 Miles (32km)	

Image EC Code: PCMM08

LENGTH 106cm (42")	HEIGHT 110cm (43")	WIDTH 63cm (25")
SEAT WIDTH 45cm (18")	MAX CAPACITY 136kg (21St)	WEIGHT 87kg (192lb)
MAXIMUM GRADIENT 8°	MAXIMUM RANGE 20 Miles (32km)	

* Dependant on user mass, terrain and battery condition.

Leaflet Code: LL148
If you would like to view this product please contact your local stockist, details across.

We have a continuous product improvement policy and consequently reserve the right to amend design and specification without prior notice. All sizes and weights stated are nominal. All errors and omissions excepted. Due to the limitations of photography and the printing process colours shown may not be 100% accurate.

Sunfire Plus EC Code: PCMM06

LENGTH 95cm (37")	HEIGHT 110cm (43")	WIDTH 60cm (24")
SEAT WIDTH 45cm (18")	MAX CAPACITY 136kg (21St)	WEIGHT 96kg (212lb)
MAXIMUM GRADIENT 8°	MAXIMUM RANGE 20 Miles (32km)	

Sunfire General Code: PCMM512

LENGTH 103cm (41")	HEIGHT 112cm (44")	WIDTH 64cm (25")
SEAT WIDTH 48cm (19")	MAX CAPACITY 180kg (28St)	WEIGHT 100kg (220lb)
MAXIMUM GRADIENT 8°	MAXIMUM RANGE 20 Miles (32km)	

Image GT Code: PCMM05

LENGTH 107cm (42")	HEIGHT 106cm (42")	WIDTH 64cm (25")
SEAT WIDTH 45cm (18")	MAX CAPACITY 160kg (25St)	WEIGHT 108kg (238lb)
MAXIMUM GRADIENT 8°	MAXIMUM RANGE 20 Miles (32km)	

POWERCHAIR RANGE

THE ULTIMATE IN DURABILITY, SAFETY & STYLE



CE



"Our excellent range of stylish new powerchairs"

For the ultimate in powerchair durability, safety and style why not experience the Drive Fixed Based Powerchair range. These powerchairs are easy to operate and extremely responsive, have a very small turning circle and are a pleasure to use.

GEO MICRO

CODE: S3500

"Compact & easy to disassemble"



- Designed for indoor use and limited outdoor use such as shopping centres
- Moulded PU flip up armrests that are width, height and angle adjustable
- Compact and manoeuvrable
- Easy to remove battery pack
- Top speed 4mph*
- On or Off board charging
- PG Drives VR2 controller
- Articulating front beam enabling easy use over gradients and uneven surfaces
- Up to 6 miles on a full battery charge*
- Weight capacity 115kg (18 stone)
- Tried and tested industry standard electronics
- Height adjustable seat and footplate
- Anti tip wheels
- Choice of red or blue

SUNFIRE PLUS GT

CODE: PCMM07

"sleek & robust design"



- Mid size rear wheel drive powerchair
- Height adjustable reclining padded seat
- Moulded PU flip up armrests that are width, height and angle adjustable
- Flip up footplate, angle and height adjustable
- PG Drives VR2 Controller
- Articulating front beam enabling easy use over gradients and uneven surfaces
- Up to 20 miles on a full battery charge*
- Sleek and robust design
- Height adjustable and removable headrest
- Powerful 350W motors
- Solid grey tyres with modern wheel design
- Weight capacity 136kg (21 stone)
- Top speed of 4mph*
- Tried and tested industry standard electronics
- Anti tip wheels

IMAGE EC

CODE: PCMM08

"Excellent manoeuvrability"



- Mid wheel drive for optimum manoeuvrability
- Six wheels for extra stability
- Height adjustable reclining padded seat
- Moulded PU flip up armrests that are width, height and angle adjustable
- Flip up footplate, angle and height adjustable
- Up to 20 miles on a full battery charge*
- Height adjustable and removable headrest
- Weight capacity 136kg (21 stone)
- Top speed of 4mph*
- Tried and tested industry standard electronics
- Anti tip wheels

SUNFIRE PLUS EC

CODE: PCMM06

"Excellent value for money"



- Mid size rear wheel drive powerchair
- Height adjustable reclining padded seat
- Moulded PU flip up armrests that are width, height and angle adjustable
- Flip up footplate, angle and height adjustable
- PG Drives VR2 Controller
- Articulating front beam enabling easy use over gradients and uneven surfaces
- Up to 20 miles on a full battery charge*
- Sleek design
- Height adjustable and removable headrest
- Solid grey tyres with modern wheel design
- Weight capacity 136kg (21 stone)
- Top speed of 4mph*
- Tried and tested industry standard electronics
- Anti tip wheels

SUNFIRE GENERAL

CODE: PCMM512

"Comfortable & accommodating large seat"



- Mid size rear wheel drive powerchair
- Height adjustable reclining padded seat
- Moulded PU flip up armrests that are width, height and angle adjustable
- Flip up footplate, angle and height adjustable
- PG Drives VR2 Controller
- Articulating front beam enabling easy use over gradients and uneven surfaces
- Innovative one-touch height adjustable armrest
- Larger seat for additional comfort
- Up to 20 miles on a full battery charge*
- Height adjustable and removable headrest
- Powerful 450W motors
- Bariatric weight capacity – 180kg (28 stone)
- Top speed of 4mph*
- Tried and tested industry standard electronics
- Anti tip wheels

IMAGE GT

CODE: PCMM05

"Improved kerb climbing ability"



- Mid wheel drive for optimum manoeuvrability
- Six wheels for extra stability
- Height adjustable reclining padded seat
- Moulded PU flip up armrests that are width, height and angle adjustable
- Flip up footplate, angle and height adjustable
- Innovative one-touch height adjustable armrest
- Powerful 450W motors
- Improved gradient climbing ability
- Compact yet sturdy design
- Height adjustable and removable headrest
- Up to 20 miles on a full battery charge*
- Weight capacity 160kg (25 stone)
- Top speed of 4mph*
- Tried and tested industry standard electronics
- Anti tip wheels