



Parts List - Spirit Wheelchair

1. Z18010-01: Centre pin assembly
2. Z18130-01*: RH cross brace assembly
3. Z18140-01*: LH cross brace assembly
4. Z18210-01: Seat canvas
5. Z18250-01: Back tension straps
6. Z18260-01: Back cushion
7. Z18180-01*: Cross brace support
8. Z17350-01: Half fold back mechanism
9. Z18325-01*: Push handle
10. Z18425-01: Castor wheel & fork ass'y
11. Z18500-01: Armpad (Left- Hand)
12. Z18550-01: Armpad (Right-Hand)
13. Z18600-01: LH legrest
14. Z18700-01: RH legrest
15. Z18725-01: LH brake
16. Z18470-01: Axle bracket
17. Z18750-01: RH brake
18. Z18450-01: Quick release pin for wheel
19. Z18460-01: 24" rear wheel
20. Z17800-01: Stepper tube c/w anti-tips
21. Z18525-01: Armrest button, spring and nut

: denotes colour suffix required