



TITLE:
LAWC001
PARTS MANUAL

DATE: 06/09/03
DRAWN: PMK
SCALE: NTS

SHEET No: 1 of 1
CHECKED: REVISION No: -
DWG No: Z16002

24	2	Z10593	ECONOMY STEEL SIDE PANEL
23	2	Z15410	CASTOR ASSEMBLY
22	2	Z10730	BRAKE CLAMP ASSEMBLY
21	4	Z10594	ARMREST SELF TAPPING SCREW
20	2	Z10592	ECONOMY ARM PAD
19	2	Z10450	QUICK RELEASE PIN ASSEMBLY
18	4	Z99965	M5 X 30 BUTTON HEAD SCREW
17	1	Z99735	M8 X 80 BUTTON HEAD HEX SKT BOLT
16	1	Z99403	M8 NYLOC NUT (T) TYPE. BZP
15	1	Z16750	TC001 BRAKE UNIT ASSEMBLY (LEFT)
14	1	Z16700	TC001 BRAKE UNIT ASSEMBLY (RIGHT)
13	1	Z16650	TC001 LHS LEG REST ASSEMBLY
12	1	Z16600	TC001 RHS LEG REST ASSEMBLY
11	2	Z16460	24" SELF PROPEL WHEEL ASSEMBLY
10	1	Z16350	HALF FOLD BACK MECHANISM (RHS)
9	1	Z16350	HALF FOLD BACK MECHANISM (LHS)
8	2	Z16325	TC001 PUSH HANDLE ASSEMBLY
7	1	Z16260	TC001 18" BACK CANVAS ASSEMBLY
6	1	Z16210	TC001 18" SEAT CANVAS ASSEMBLY
5	1	Z16140	TC001 LHS CROSS BRACE ASSEMBLY
4	1	Z16130	TC001 RHS CROSS BRACE ASSEMBLY
3	1	Z16035	LAWC001 LHS FRAME ASSEMBLY
2	1	Z16025	LAWC001 RHS FRAME ASSEMBLY
1	1	Z16010	CENTRE PIN ASSEMBLY
Item	Qty	Name	Description
Parts List			



TITLE:

LAWC001 WHEELCHAIR
PARTS MANUAL

DATE:

06/09/03

SHEET No:

2 OF 2

DRAWN:

PMK

CHECKED:

REVISION No:

—

SCALE:

NTS

DWG No:

Z16002