

MEDICARE TECHNOLOGY: WHEELCHAIR / SCOOTER RECOMMENDED SERVICE ITEMS

	Enigma Steel	Travel Chair	Enigma XS	Energi Powerchair	Mambo Powerchair	Travelease	Pioneer	Rio	Monte Carlo	Apollo	M36 / M44 / M46	M48		Enigma Steel	Travel Chair	Enigma XS	Energi Powerchair	Mambo Powerchair	Travelease	Pioneer	Rio	Monte Carlo	Apollo	M36 / M44 / M46	M48																		
Section 1: Controller														Section 5: Motor/s																													
On/Off Switch																											Wiring																
Joystick / Wigwag																											Mounting																
Joystick Gaitor																											Function																
Output Plug																											Noise Level																
Seal																											Commutator Cleaned																
Operation																											Brushes																
Dynamic Braking																											Section 6: Upholstery																
Static Braking																											Seat																
Section 2: Batteries														Section 7: Body																													
Type (SLA/AGM/Gel)																											Mounting																
Compartment																											Condition																
Connections																											Seat Mechanism																
Wirings																											Seat Adjustment																
Battery Test																											Section 8: Electrics																
Section 3: Wheels														Section 9: Drive Train																													
Tyre Pressure																											Transaxle Seals																
Tyre Wear																											Transaxle Noise Lever																
Wheel Bearings																											Transaxle Mount Brackets																
Hub Caps / Disks																											Section 10: Test Run																
Wheel Nuts Tight																											Test Run																
Tab Washers																											Section 10: Test Run																
Lubricated / Greased																											Test Run																
Section 4: Chassis / Frame														Section 10: Test Run																													
General Condition																											Test Run																
Castor Assemblies																											Section 10: Test Run																
Manual Brakes																											Test Run																
Suspension																											Section 10: Test Run																
Steering Linkage																											Test Run																
Swivel Bearings																											Section 10: Test Run																
Tiller Pivot																											Test Run																
Handbrake																											Section 10: Test Run																
Footrest Lock Mechanism																											Test Run																