



TITLE:
**ENIGMA TRAVEL CHAIR
 PARTS MANUAL**

DATE: 06/09/03
 DRAWN: PMK
 SCALE: NTS

SHEET No: 1 of 1
 CHECKED:
 REVISION No: -
 DWG No: Z16001

22	1	Z99403	M8 NYLOC NUT (T) TYPE. BZP
21	1	Z99735	M8 X 80 BUTTON HEAD HEX SKT BOLT
20	2	Z10730	BRAKE CLAMP ASSEMBLY
19	1	Z16750	TC001 BRAKE UNIT ASSEMBLY (LEFT)
18	1	Z16700	TC001 BRAKE UNIT ASSEMBLY (RIGHT)
16	4	Z99965	M5 X 30 BUTTON HEAD SCREW
15	2	Z10592	ECONOMY ARM PAD
14	1	Z16010	CENTRE PIN ASSEMBLY
13	2	Z15410	CASTOR ASSEMBLY
12	1	Z16650	TC001 LHS LEG REST ASSEMBLY
11	1	Z16600	TC001 RHS LEG REST ASSEMBLY
10	2	Z16480	TC001 REAR WHEEL ASSEMBLY
9	1	Z16350	HALF FOLD BACK MECHANISM (RHS)
8	1	Z16350	HALF FOLD BACK MECHANISM (LHS)
7	2	Z16325	TC001 PUSH HANDLE ASSEMBLY
6	1	Z16260	TC001 18" BACK CANVAS ASSEMBLY
5	1	Z16210	TC001 18" SEAT CANVAS ASSEMBLY
4	1	Z16140	TC001 LHS CROSS BRACE ASSEMBLY
3	1	Z16130	TC001 RHS CROSS BRACE ASSEMBLY
2	1	Z16030	TC001 LHS FRAME ASSEMBLY
1	1	Z16020	TC001 RHS FRAME ASSEMBLY
Item	Qty	Name	Description
Z16001			



TITLE:
ENIGMA TRAVEL CHAIR
PARTS MANUAL

DATE: 06/09/03
DRAWN: PMK
SCALE: NTS

SHEET No: 2 OF 2
REVISION No: -
DWG No: Z16001