



TITLE:  
LAWC002  
PARTS MANUAL

DATE: 10/11/05  
DRAWN: PMK  
SCALE: NTS

SHEET No: 1 of 1  
CHECKED: REVISION No: -  
DWG No: Z16003

|                   |            |             |                                  |
|-------------------|------------|-------------|----------------------------------|
| 22                | 2          | Z10730      | 7/8" BRAKE CLAMPS (NOT SHOWN)    |
| 23                | 1          | EB003       | DUAL EASIBRAKE (NOT SHOWN)       |
| 1                 | 1          | Z16010      | CENTRE PIN ASSEMBLY              |
| 2                 | 1          | Z16040      | LAWC002 RHS FRAME ASSEMBLY       |
| 3                 | 1          | Z16045      | LAWC002 LHS FRAME ASSEMBLY       |
| 4                 | 1          | Z16130      | TC001 RHS CROSS BRACE ASSEMBLY   |
| 5                 | 1          | Z16140      | TC001 LHS CROSS BRACE ASSEMBLY   |
| 6                 | 1          | Z16210      | TC001 18" SEAT CANVAS ASSEMBLY   |
| 7                 | 1          | Z16260      | TC001 18" BACK CANVAS ASSEMBLY   |
| 8                 | 2          | Z16330      | EXTENDED PUSH HANDLE ASSEMBLY    |
| 9                 | 1          | Z16350      | HALF FOLD BACK MECHANISM (LHS)   |
| 10                | 1          | Z16350      | HALF FOLD BACK MECHANISM (RHS)   |
| 11                | 2          | Z10480      | 12-1/2" TRANSIT WHEEL ASSEMBLY   |
| 12                | 1          | Z16600      | TC001 RHS LEG REST ASSEMBLY      |
| 13                | 1          | Z16650      | TC001 LHS LEG REST ASSEMBLY      |
| 14                | 2          | Z16450      | QUICK RELEASE PIN ASSEMBLY       |
| 16                | 1          | Z99403      | M8 NYLOC NUT (T) TYPE. BZP       |
| 17                | 1          | Z99735      | M8 X 80 BUTTON HEAD HEX SKT BOLT |
| 18                | 4          | Z99965      | M5 X 30 BUTTON HEAD SCREW        |
| 15                | 2          | Z10592      | ECONOMY ARM PAD                  |
| 19                | 4          | Z10594      | ARMREST SELF TAPPING SCREW       |
| 20                | 2          | Z15410      | CASTOR ASSEMBLY                  |
| 21                | 2          | Z10593      | ECONOMY STEEL SIDE PANEL         |
| <b>Item</b>       | <b>Qty</b> | <b>Name</b> | <b>Description</b>               |
| <b>Parts List</b> |            |             |                                  |



TITLE:

LAWC002 WHEELCHAIR  
PARTS MANUAL

DATE:

10/11/05

SHEET No:

2 OF 2

DRAWN:

PMK

CHECKED:

REVISION No:

—

SCALE:

NTS

DWG No:

Z16003