



TC002 Parts List

Drive Medical Ltd,
Ainley's Industrial Estate,
Elland, West Yorkshire HX5 9JP.
Telephone: +44 (0) 1422 314488
Email: technical@drivemedical.co.uk
Web: www.drivemedical.co.uk/technical

Document Code: Z16004

No.	Qty	Code	Description
1.	1	Z16050-01	TC002 RHS frame assembly
2.	1	Z16055-01	TC002 LHS frame assembly
3.	1	Z16150-01	TC002 (19") RHS cross brace assembly
4.	1	Z16160-01	TC002 (19") LHS cross brace assembly
5.	2	Z16180-01	Cross brace support bracket
6.	1	Z16215-01	TC002 (19") seat canvas assembly
7.	1	716	Lap belt
8.	1	Z16265-01	TC002 (19") back canvas assembly
9.	1	Z16340-01	TC002 RHS push handle
10.	1	Z16345-01	TC002 LHS push handle
11.	1	Z16360-01	RHS half-fold back mechanism
12.	1	Z16350-01	LHS half-fold back mechanism
13.	2	Z16400-01	6" front castor assembly
14.	2	Z16490-01	5-spoke rear wheel assembly
15.	1	Z16600-01A	TC002 RHS leg rest assembly
16.	1	Z16650-01A	TC002 LHS leg rest assembly
17.	1	Z16705-01	RHS fixed brake
18.	1	Z16755-01	LHS fixed brake
19.	2	Z16800-01	TC002 stepper tube
20.	2	Z10531-01NY	Nylon upholstered arm pad